

HU-FRIEDY® KEY OPINION LEADER

DR. BORA KORKUT

RESTORATIVE & ISOLATION KIT



Dr. Bora Korkut,
D.D.S., Ph.D.

Academic

2008 Dental Licence [DDS] (Marmara University, Faculty of Dentistry)

2015 Doctorate Degree [PhD] (Marmara University, Institute of Health Sciences)

2015 Assistant Professor / Restorative Dentistry Specialist (Marmara University, Department of Restorative Dentistry)

Fields of Interest

Aesthetic & Restorative Dentistry

Digital Dentistry

Direct Resin Restorations

Cad/Cam Restorations

Macro Dental Photography

Professional Memberships

Turkish Dental Association (Member)

Turkish Society of Restorative Dentistry (Member)

Turkish Academy of Esthetic Dentistry (Member)

Digital Dentistry Society (Member)

Megagen International Network of Education

and Clinical Research (Member)

Mail: bora.korkut@marmara.edu.tr

Instagram: dr.borakorkut



BORA KORKUT
ESTHETIC & RESTORATIVE DENTISTRY



It's not only the skill of the doctor but also the type of the instrument that determines the quality of restoration. During the restoration, it is impossible for me to simulate anatomical details precisely without these specially designed instruments.

Dr. Bora Korkut



RESTORATIVE KIT


This restoration kit offers excellent manipulation control for both anterior and posterior direct resin restorations.

IME14DIN5_
IMS™ Infinity Series™
Cassette DIN, 5 Instruments

CHOOSE THE COLOR

Just add the color coding number to the last part of the partcode (replace “_” with the number)

Red	●	0
Grey	●	1
White	●	2
Orange	●	3
Purple	●	4
Yellow	●	5
Lavendar	●	6
Black	●	7
Blue	●	8
Green	●	9
Ocean Blue	●	11



WHY YOU SHOULD USE AN IMS™ INFINITY SERIES™ CASSETTE SYSTEM:

By standardising your instruments with IMS™ System, additional time and revenue is generated while extending the useful life of your Hu-Friedy® instruments. Dental professionals state that by implementing IMS™ System, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.

*When compared to single instrument decontamination. Based on market survey results.

Visit us online at HuFriedy.eu

Hu-Friedy is now a proud member of

HuFriedyGroup
The Best In Practice

RESTORATIVE KIT

This restoration kit offers excellent manipulation control for both anterior and posterior direct resin restorations.




CVKOT1X
XTS™ PFI Essential
WeRestore.it Spatula, grey



TNFF1/2
Ferraris Main shaper



AKRO-FLEX
Akro-Flex™ Composite
Instrument



TNPCCI
Plastic Filling Instrument Duck
Head hdl #6S XTS




TNFF5/6
Ferraris 3D Shaper

RUBBER DAM KIT


Rubber-Dam isolation is a 'must' for direct anterior functional and aesthetic restorations. And Hu-Friedy Black Line clamps are perfect tools for optimum isolation, stabilization and gingival retraction. These black clamps are also unique for macro dental photography and digital scanning as to prevent unwilling flash light reflections.

BLACK LINE RUBBER DAM CLAMPS


- Satin Steel® which resists corrosion while aiding in flexibility and memory
- A dark coating allowing for improved intraoral contrast and visual acuity
- Matte finish for decreased light reflection and glare




RDCM00X (2x)
Rubber Dam Clamps
#00 UK anterior,
Black Line




RDCM2X (2x)
Rubber Dam Clamps
#2 UK premolar,
Black Line



RDCM212X
Rubber Dam Clamps
#212 anterior,
Black Line



RDCM212SAX
Rubber Dam Clamps
#212SA anterior,
Black Line



RDCM9X
Rubber Dam Clamps
#9 anterior with wings,
Black Line




RDAF6
Rubber Dam Frame, Adult



IMS-1273
IMS™ Accessory Box, large



RDF
Rubber Dam Forceps



RDP
Rubber Dam Punch

CHOOSE THE COLOR

Just add the color coding number to the last part of the partcode (replace “_” with the number)

Red	●	0
Grey	●	1
White	●	2
Orange	●	3
Purple	●	4
Yellow	●	5
Lavendar	●	6
Black	●	7
Blue	●	8
Green	●	9
Ocean Blue	●	11

IMEDIN16_

IMS™ Infinity Series™ Cassette DIN,
16 Instruments



ADDITIONALS

This restoration kit offers excellent manipulation control for both anterior and posterior direct resin restorations.



PROGS

Chu Gauges
Proportion Kit
hdl #6S
(1 PROCHDLS,
2 T-BAR +
2 INLINE tips)



INLINEREF

Chu Gauges INLINE Refill Tips
3pcs./pack.



TBARREF

Chu Gauges T-BAR refill tips
3pcs./pack.



TNCIGFTMI4

Plastic Filling Instrument Goldstein #4 hdl
#6S XTS™, Flexi-Thin mini



TNCALIBRA

XTS™ PFI WeRestore.it
Calibra™ Instrument



10-130-05E

Scalpel Handle #5E, European handle



MIR5HDBE

High Definition Mirror #5
single-sided, European,
Black Line, 1pc./pack



MHE6B

Mirror Handle hdl #6
European Thread,
Black Line



S17SCX

Scissors Iris #17 straight SuperCut
11,5cm Black Line

Free Call: 00800 48 37 43 39 | Free Fax: 00800 48 37 43 40

E-Mail: info@hu-friedy.eu | Website: HuFriedy.eu

Facebook: @HuFriedyEU | Instagram: @hu_friedy_europe

LinkedIn: @Hu-Friedy Mfg. Co., LLC.EU

All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC.,
its affiliates or related companies, unless otherwise noted.
©2021 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-381GB/0421

Hu-Friedy is now a proud member of

HuFriedyGroup
The Best In Practice